



RPL23 Antibody

Product Code	CSB-PA070295
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62829
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein L23.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPL23; 60S ribosomal protein L23; 60S ribosomal protein L17
Immunogen Species	Homo sapiens (Human)
Target Names	RPL23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.