



SLC22A6 Antibody

Product Code	CSB-PA070296
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95742
Immunogen	Synthesized peptide derived from the Internal region of Human OAT1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	Solute carrier family 22 member 6; Organic anion transporter 1; Hoat1; Renal organic anion transporter 1; hROAT1; PAH transporter; hPAHT; SLC22A6; OAT1; PAHT
Immunogen Species	Homo sapiens (Human)
Target Names	SLC22A6