



Phospho-CD28 (Y218) Antibody

Product Code	CSB-PA070302
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10747
Immunogen	Synthesized peptide derived from Human CD28 around the phosphorylation site of Y218.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD28; T-cell-specific surface glycoprotein CD28; TP44; CD antigen CD28
Immunogen Species	Homo sapiens (Human)
Target Names	CD28
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.