




Phospho-CD40 (T254) Antibody

Product Code	CSB-PA080031
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25942
Immunogen	Synthesized peptide derived from Human CD40 around the phosphorylation site of T254.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD40; TNFRSF5; Tumor necrosis factor receptor superfamily member 5; B-cell surface antigen CD40; Bp50; CD40L receptor; CDw40; CD antigen CD40
Immunogen Species	Homo sapiens (Human)
Target Names	CD40

Image



Western Blot analysis of COLO cells using Phospho-CD40 (T254) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.