



MAD2L1BP Antibody

Product Code	CSB-PA080033
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15013
Immunogen	Synthesized peptide derived from the N-terminal region of Human MAD2L1BP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAD2L1BP; CMT2; KIAA0110; MAD2L1-binding protein; Caught by MAD2 protein
Immunogen Species	Homo sapiens (Human)
Target Names	MAD2L1BP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.