



CELSR1 Antibody

Product Code	CSB-PA080041
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ6
Immunogen	Synthesized peptide derived from the Internal region of Human CDHF9.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CELSR1; CDHF9; FMI2; Cadherin EGF LAG seven-pass G-type receptor 1; Cadherin family member 9; Flamingo homolog 2; hFmi2
Immunogen Species	Homo sapiens (Human)
Target Names	CELSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.