



SH3PXD2A Antibody

Product Code	CSB-PA080044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCZ1
Immunogen	Synthesized peptide derived from the N-terminal region of Human Fish.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SH3PXD2A; FISH; KIAA0418; SH3MD1; TKS5; SH3 and PX domain-containing protein 2A; Adapter protein TKS5; Five SH3 domain-containing protein; SH3 multiple domains protein 1; Tyrosine kinase substrate with five SH3 domains
Immunogen Species	Homo sapiens (Human)
Target Names	SH3PXD2A