



Phospho-MDC1 (S513) Antibody

Product Code	CSB-PA080045
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14676
Immunogen	Synthesized peptide derived from Human MDC1 around the phosphorylation site of S513.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MDC1; KIAA0170; NFBD1; Mediator of DNA damage checkpoint protein 1; Nuclear factor with BRCT domains 1
Immunogen Species	Homo sapiens (Human)
Target Names	MDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.