



RAB11FIP3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA080054 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75154 |
| Immunogen | Synthesized peptide derived from the C-terminal region of Human Rab11-FIP3. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | RAB11FIP3; ARFO1; KIAA0665; Rab11 family-interacting protein 3; FIP3-Rab11; Rab11-FIP3; Arfophilin-1; EF hands-containing Rab-interacting protein; Eferin; MU-MB-17.148 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RAB11FIP3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |