



Phospho-ACAP1 (S554) Antibody

Product Code	CSB-PA080059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15027
Immunogen	Synthesized peptide derived from Human ACAP1 around the phosphorylation site of S554.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACAP1; CENTB1; KIAA0050; Arf-GAP with coiled-coil; ANK repeat and PH domain-containing protein 1; Centaurin-beta-1; Cnt-b1
Immunogen Species	Homo sapiens (Human)
Target Names	ACAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.