



TLK1 Antibody

Product Code	CSB-PA080080
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4F6
Immunogen	Synthesized peptide derived from Human TLK1 around the non-phosphorylation site of S764.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TLK1; KIAA0137; Serine/threonine-protein kinase tousled-like 1; PKU-beta; Tousled-like kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	TLK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.