



Phospho-CDC6 (S54) Antibody

Product Code	CSB-PA080087
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99741
Immunogen	Synthesized peptide derived from Human Cdc6 around the phosphorylation site of S54.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDC6; CDC18L; Cell division control protein 6 homolog; CDC6-related protein; Cdc18-related protein; HsCdc18; p62(cdc6); HsCDC6
Immunogen Species	Homo sapiens (Human)
Target Names	CDC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.