



Phospho-CDC25A (S178) Antibody

Product Code	CSB-PA080094
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30304
Immunogen	Synthesized peptide derived from Human Cdc25A around the phosphorylation site of S178.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDC25A; M-phase inducer phosphatase 1; Dual specificity phosphatase Cdc25A
Immunogen Species	Homo sapiens (Human)
Target Names	CDC25A