



CDC25C Antibody

Product Code	CSB-PA080109
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30307
Immunogen	Synthesized peptide derived from Human Cdc25C around the non-phosphorylation site of T48.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDC25C; M-phase inducer phosphatase 3; Dual specificity phosphatase Cdc25C
Immunogen Species	Homo sapiens (Human)
Target Names	CDC25C