




Mono-Methyl-Histone H3 (Lys9) Antibody

Product Code	CSB-PA080120
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431/Q71DI3/P84243
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)

Image



Western blot analysis of 1) HeLa, 2) Raw264.7, 3) Rat Testis Tissue, diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.