




Tri-Methyl-Histone H3 (Lys14) Antibody

Product Code	CSB-PA080124
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431/Q71DI3/P84243
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H3A/HIST2H3A/H3F3A

Image



Western blot analysis of 1) HeLa, 2) Rat Testis Tissue, diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.