




# ATP1B1 Antibody

<b>Product Code</b>	CSB-PA080138
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P05026
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse,Rat,Human
<b>Tested Applications</b>	WB;WB:1:1000-2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	Sodium/potassium-transporting ATPase subunit beta-1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ATP1B1

## Image



Western blot analysis of 1) Mouse Brain, 2) Rat Brain, diluted at 1:2000.