




CTNNB1 Antibody

Product Code	CSB-PA080139
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB;WB:1:2000-4000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNB1

Image



Western blot analysis of 1) 3T3, 2) Rat Liver, diluted at 1:2,000.