




CTNNB1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA080139 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P35222 |
| Immunogen | Recombinant Protein |
| Raised In | Rabbit |
| Species Reactivity | Mouse,Rat |
| Tested Applications | WB;WB:1:2000-4000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | CTNNB1 |

Image



Western blot analysis of 1) 3T3, 2) Rat Liver, diluted at 1:2,000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.