




PEP1 Antibody

Product Code	CSB-PA080140
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04711
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Maize
Tested Applications	WB;WB:1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Maize
Target Names	PEP1

Image



Western blot analysis of Maize, diluted at 1:4000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.