




# PGR Monoclonal Antibody

<b>Product Code</b>	CSB-MA080142
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P06401
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,IHC,IF;Recommended dilution:IHC:1:100-1:200,IF:1:100-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	Progesterone receptor
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PGR
<b>Clone No.</b>	3G5

## Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma using PR Monoclonal Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.