





ZBTB45 Antibody

Product Code	CSB-PA080144
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96K62
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Zinc finger and BTB domain-containing protein 45
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB45

Image



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.