





COX4I1 Antibody

Product Code	CSB-PA080145
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13073
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000, IHC:1:500-1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Immunogen Species	Homo sapiens (Human)
Target Names	COX4I1

Image



Immunohistochemical analysis of paraffin-embedded human Colon carcinoma using COX4 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) Mouse Brain, 4) Rat Brain using COX4 Polyclonal Antibody.