





CDH1 Antibody

Product Code	CSB-PA080147
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12830
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC; WB:1:1000, IHC:1:500-1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Cadherin-1
Immunogen Species	Homo sapiens (Human)
Target Names	CDH1

Image



Immunohistochemical analysis of paraffin-embedded human Lung carcinoma using E-cadherin Polyclonal Antibody.



Western blot analysis of HeLa using E-cadherin Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.