





HSP90AA1 Antibody

Product Code	CSB-PA080148
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07900
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Heat shock protein HSP 90-alpha
Immunogen Species	Homo sapiens (Human)
Target Names	HSP90AA1

Image



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using HSP90A Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) Mouse Brain, 4) Rat Brain using HSP90A Polyclonal Antibody.