





YWHAE Antibody

Product Code	CSB-PA080149
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62258
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	14-3-3 protein epsilon
Immunogen Species	Homo sapiens (Human)
Target Names	YWHAE

Image



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using 14-3-3 ε Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain, 2) Rat Brain using 14-3-3 ε Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.