





ARNT Antibody

Product Code	CSB-PA080154
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27540
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Aryl hydrocarbon receptor nuclear translocator
Immunogen Species	Homo sapiens (Human)
Target Names	ARNT

Image



Immunohistochemical analysis of paraffin-embedded human Tonsil Tissue using HIF-1 β Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain, 2) Rat Brain using HIF-1 β Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.