





MAPK7 Antibody

Product Code	CSB-PA080163
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13164
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Mitogen-activated protein kinase 7
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK7

Image



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using ERK 5 Polyclonal Antibody.



Western blot analysis of 1) Hela, 2) 293T, 3) Mouse Skeletal Muscle, 4) Rat Kidney, 5) Rat Skeletal Muscle using ERK 5 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.