





SNAI2 Antibody

Product Code	CSB-PA080165
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43623
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Zinc finger protein SNAI2
Immunogen Species	Homo sapiens (Human)
Target Names	SNAI2

Image



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using SLUG Polyclonal Antibody.



Western blot analysis of 1) MCF7, 2) Mouse Heart Tissue, 3) Rat Heart Tissue, 4) Rat Brain Tissue using SLUG Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.