





GSK3B Antibody

Product Code	CSB-PA080166
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49841
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Glycogen synthase kinase-3 beta
Immunogen Species	Homo sapiens (Human)
Target Names	GSK3B

Image



Immunohistochemical analysis of paraffin-embedded Mouse Lung Tissue using GSK3 β Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Mouse Heart Tissue, 3) Rat Brain Tissue, 4) Rat Heart Tissue using GSK3 β Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.