





# GSK3B Antibody

<b>Product Code</b>	CSB-PA080166
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P49841
<b>Immunogen</b>	Recombinant Protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse,Rat
<b>Tested Applications</b>	WB, IHC; WB:1:1000-2000, IHC:1:200-500
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	Glycogen synthase kinase-3 beta
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GSK3B

## Image



Immunohistochemical analysis of paraffin-embedded Mouse Lung Tissue using GSK3β Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Mouse Heart Tissue, 3) Rat Brain Tissue, 4) Rat Heart Tissue using GSK3β Polyclonal Antibody.