




ECFP-Tag Monoclonal Antibody

Product Code	CSB-MA080170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	All
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0
Product Type	Monoclonal Antibody
Immunogen Species	All
Target Names	ECFP-Tag
Clone No.	18F5

Image



Western blot analysis of Recombinant ECFP Protein using ECFP-Tag Monoclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.