





TP53 Antibody

Product Code	CSB-PA080174
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:500-1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Cellular tumor antigen p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Immunohistochemical analysis of paraffin-embedded Mouse Colon Tissue using p53 Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue using p53 Polyclonal Antibody.