





TP53 Antibody

Product Code	CSB-PA080174
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:500-1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Cellular tumor antigen p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Immunohistochemical analysis of paraffin-embedded Mouse Colon Tissue using p53 Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue using p53 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.