





Histone H2A.X Antibody

Product Code	CSB-PA080179
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16104
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Histone H2AX
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFX

Image



Immunohistochemical analysis of paraffin-embedded Mouse Colon Tissue using Histone H2A.X Polyclonal Antibody.



Western blot analysis of 1) 3T3, 2) Rat Brain Tissue using Histone H2A.X Polyclonal Antibody.