





Histone H2A.X Antibody

Product Code	CSB-PA080179
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16104
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Histone H2AX
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFX

Image



Immunohistochemical analysis of paraffin-embedded Mouse Colon Tissue using Histone H2A.X Polyclonal Antibody.



Western blot analysis of 1) 3T3, 2) Rat Brain Tissue using Histone H2A.X Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.