




BIRC5 Antibody

Product Code	CSB-PA080180
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15392
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:200-500, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Baculoviral IAP repeat-containing protein 5
Immunogen Species	Homo sapiens (Human)
Target Names	BIRC5

Image



Immunohistochemical analysis of paraffin-embedded human Colon carcinoma using Survivin Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.