





NMBR Antibody

Product Code	CSB-PA080188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28336
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Neuromedin-B receptor
Immunogen Species	Homo sapiens (Human)
Target Names	NMBR

Image



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using NMBR Polyclonal Antibody.



Western blot analysis of 1) Human Brain Tissue, 2) Rat Brain Tissue using NMBR Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.