





TF Antibody

Product Code	CSB-PA080189
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02787
Immunogen	Purified Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:5000-10000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Serotransferrin
Immunogen Species	Homo sapiens (Human)
Target Names	TF

Image



Immunohistochemical analysis of paraffin-embedded Human Liver Tissue using TF Polyclonal Antibody.



Western blot analysis of Human Blood sample using TF Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.