





CASR Antibody

Product Code	CSB-PA080191
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41180
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Extracellular calcium-sensing receptor
Immunogen Species	Homo sapiens (Human)
Target Names	CASR

Image



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using CASR Polyclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue using CASR Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.