





CTNNB1 Monoclonal Antibody

Product Code	CSB-MA080197
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Catenin beta-1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNB1
Clone No.	12E8

Image



Immunohistochemical analysis of paraffin-embedded Human Lung carcinoma using Catenin-β Monoclonal Antibody.



Western blot analysis of 1) HeLa, 2) 293T, 3) Mouse Liver Tissue, 4) Rat Liver Tissue using Catenin-β Monoclonal Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.