




PIK3R1 Antibody

Product Code	CSB-PA080202
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27986
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Rat,Mouse
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Phosphatidylinositol 3-kinase regulatory subunit alpha
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R1

Image



Western blot analysis of 1) 3T3, 2) Raw264.7, 3) Rat Kidney Tissue, 4) Rat Liver Tissue using PI 3-kinase p85 α Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.