



LexA DNA Binding Region Antibody

Product Code	CSB-PA080206
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Tested Applications	WB;WB:1:10000-50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0

Image



Western blot analysis of Recombinant LexA protein using LexA DNA Binding Region Polyclonal Antibody.