





Histone H3 Antibody

Product Code	CSB-PA080208
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Histone H3.1
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H3A

Image



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using Histone H3 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3 using Histone H3 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.