





Histone H3 Antibody

Product Code	CSB-PA080208
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:500-1000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Histone H3.1
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H3A

Image



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using Histone H3 Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3 using Histone H3 Polyclonal Antibody.