





PIK3R1 Antibody

Product Code	CSB-PA080209
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27986
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Phosphatidylinositol 3-kinase regulatory subunit alpha
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R1

Image



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using PI 3-kinase p85α Polyclonal Antibody.



Western blot analysis of 1) HepG2, 2) 293T using PI 3-kinase p85α Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.