





STAT2 Antibody

Product Code	CSB-PA080210
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52630
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Signal transducer and activator of transcription 2
Immunogen Species	Homo sapiens (Human)
Target Names	STAT2

Image



Immunohistochemical analysis of paraffin-embedded Human Lung carcinoma using Stat2 Polyclonal Antibody.



Western blot analysis of 1) K562, 2) Mouse Heart Tissue, 3) Rat Heart Tissue using Stat2 Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.