





BIRC5 Monoclonal Antibody

Product Code	CSB-MA080214
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15392
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Baculoviral IAP repeat-containing protein 5
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BIRC5
Clone No.	23F1

Image



Immunohistochemical analysis of paraffin-embedded human colon carcinoma using Survivin Monoclonal Antibody.



Western blot analysis of 1) Hela, 2) 293, 3) PC12 using Survivin Monoclonal Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.