





# PIK3R2 Antibody

<b>Product Code</b>	CSB-PA080215
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00459
<b>Immunogen</b>	Recombinant Protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, IHC; WB:1:1000-2000, IHC:1:50-200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	Phosphatidylinositol 3-kinase regulatory subunit beta
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PIK3R2

## Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma using PI 3-Kinase p85 $\beta$  Polyclonal Antibody.



Western blot analysis of HepG2 using PI 3-Kinase p85 $\beta$  Polyclonal Antibody.