





CANX Antibody

Product Code	CSB-PA080219
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27824
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC; WB:1:1000-2000, IHC:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	Calnexin
Immunogen Species	Homo sapiens (Human)
Target Names	CANX

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Calnexin Polyclonal Antibody.



Western blot analysis of 1) HeLa, 2) 3T3, 3) PC12 using Calnexin Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.